

产品名称: Rab L5 Rabbit Polyclonal Antibody
产品货号: APRab16765



产品概述 (Summary)

产品名称 (Production Name)	Rab L5 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RABL5
别名 (Alternative Names)	RABL5; Rab-like protein 5
基因 ID (Gene ID)	64792.0
蛋白 ID (SwissProt ID)	Q9H7X7. The antiserum was produced against synthesized peptide derived from human RABL5. AA range: 128-177

产品应用 (Application)


稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)



研究领域 (Research Area)

图片 (Image Data)



注意事项 (Note)

Web : <https://www.enkilife.cn>

E-mail : order@enkilife.cn (销售) tech@enkilife.cn (技术支持)

Tel : 027-87002838