

**产品名称:** RAB13 Rabbit Polyclonal Antibody

**产品货号:** APRab16773

### 产品概述 (Summary)

|                   |                                                                                          |
|-------------------|------------------------------------------------------------------------------------------|
| 描述(Description)   | Rabbit polyclonal Antibody                                                               |
| 宿主(Host)          | Rabbit                                                                                   |
| 应用(Application)   | WB,ELISA                                                                                 |
| 种属反应(Reactivity)  | Human,Rat,Mouse                                                                          |
| 偶联物(Conjugation)  | Unconjugated                                                                             |
| 修饰(Modification)  | Unmodified                                                                               |
| 同种型(Isotype)      | IgG                                                                                      |
| 克隆(Clonality)     | Polyclonal                                                                               |
| 剂型(Form)          | Liquid                                                                                   |
| 保存条件(Storage)     | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| 储存溶液(Buffer)      | Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.                |
| 纯化 (Purification) | Affinity purification                                                                    |

### 产品应用(Application)

|                       |                                      |
|-----------------------|--------------------------------------|
| 稀释比(Dilution Ratio)   | WB 1:500-1:2000,ELISA 1:5000-1:20000 |
| 分子量(Molecular Weight) | 22kDa                                |

### 抗原信息(Antigen Information)

|                       |                                                                      |
|-----------------------|----------------------------------------------------------------------|
| 基因名(Gene Name)        | RAB13 GIG4                                                           |
| 别名(Alternative Names) |                                                                      |
| 基因 ID(Gene ID)        | 5872.0                                                               |
| SwissProt ID          | P51153                                                               |
| 免疫原(Immunogen)        | Synthesized peptide derived from human protein . at AA range: 40-120 |

### 研究背景 (Background)

This gene is a member of the Rab family of small G proteins and plays a role in regulating membrane trafficking between trans-Golgi network (TGN) and recycling endosomes (RE). The encoded protein is involved in the assembly of tight junctions, which are components of the apical junctional complex (AJC) of epithelial cells. The AJC plays a role in forming a barrier between luminal contents and the underlying tissue. Additional functions associated with the protein include endocytic recycling of

occludin, regulation of epithelial cell scattering, neuronal regeneration and regulation of neurite outgrowth. Alternately spliced transcript variants have been observed for this gene. A pseudogene associated with this gene is located on chromosome 12. [provided by RefSeq, Jan 2013],function:Could participate in polarized transport, in the assembly and/or the activity of tight junctions.,similarity:Belongs to the small GTPase superfamily. Rab family.,subcellular location:Cytoplasmic tight junctions or associated with vesicles scattered throughout the cytoplasm in cells lacking tight junctions.,tissue specificity:Detected in several types of epithelia, including intestine, kidney, liver and in endothelial cells.,


### 研究领域 (Research Area)

Tight junction;

### 注意事项 (Note)

For research use only.

### 图片 (Image Data)



Western blot analysis of lysates from DU145 cells, RAB13 Rabbit Polyclonal Antibody was diluted at 1:1000, 4°over night

## EnkiLife 优势产品推荐

| 产品类别                        | 产品货号     | 产品名称                                               |
|-----------------------------|----------|----------------------------------------------------|
| WB 解决方案                     | RA10020  | 2h 极速 WB 即用型全流程试剂盒                                 |
|                             | RA10021  | 4h 快速 WB 即用型全流程试剂盒                                 |
|                             | RA10042  | 5h 畅享版 WB 全流程试剂盒                                   |
|                             | RA10037  | 校准级彩色预染蛋白 Marker (8-180kDa)                        |
|                             | RA10038  | 校准级彩色预染蛋白 Marker (10-250kDa)                       |
|                             | RA10039  | 校准级高分子彩色预染蛋白 Marker (25-400kDa)                    |
| TSA 多重荧光染色试剂盒               | RA10008  | TSA 双标三色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10009  | TSA 三标四色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10010  | TSA 四标五色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10011  | TSA 五标六色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10012  | TSA 六标七色多重荧光染色试剂盒 (mIHC)                           |
| IHC 检测试剂盒                   | RA10006  | HRP Anti-Mouse/Rabbit IHC Detection System         |
|                             | RA10007  | Polymer-HRP Anti-Mouse/Rabbit IHC Detection System |
| 抗体标记试剂盒                     | RE80004p | 辣根过氧化物酶(HRP)抗体标记试剂盒                                |
|                             | RE80002q | Sulfo-NHS-生物素标记试剂盒                                 |
|                             | RE80007p | Cy3 荧光素标记试剂盒                                       |
|                             | RE80011p | Fluor488 荧光素标记试剂盒                                  |
|                             | RE80017p | Fluor750 荧光素标记试剂盒                                  |
|                             | RE80005p | 藻红蛋白(R-PE) 抗体快速标记试剂盒                               |
|                             | RE80040  | PE-Cy7 串联染料抗体快速标记试剂盒                               |
| 稳转细胞系构建服务<br>(免费赠送全膜 WB 验证) | TS-0001  | 过表达稳转细胞系构建                                         |
|                             | TS-0002  | 敲低稳转细胞系构建                                          |
|                             | TS-0003  | 敲除细胞系构建                                            |

限于篇幅，更多产品详见官网：<https://www.enkilife.cn>