

产品名称: RAB28 Rabbit Polyclonal Antibody
产品货号: APRab16783



产品概述 (Summary)

产品名称 (Production Name)	RAB28 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAB28
别名 (Alternative Names)	
基因 ID (Gene ID)	9364.0
蛋白 ID (SwissProt ID)	P51157.Synthesized peptide derived from human protein . at AA range: 80-160

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	24kDa

研究背景 (Background)

This gene encodes a member of the Rab subfamily of Ras-related small GTPases. The encoded protein may be involved in regulating intracellular trafficking. Alternative splicing results in multiple transcript variants. Pseudogenes of this gene are


产品名称: RAB28 Rabbit Polyclonal Antibody
产品货号: AP Rab16783



found on chromosomes 9 and X. [provided by RefSeq, Apr 2009],similarity:Belongs to the small GTPase superfamily. Rab family.,tissue specificity:Isoform S is detected in most tissues investigated (cortex, liver, kidney, skeletal muscle, adipose tissue, testis and urothelium), whereas isoform L is predominant in testis.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4° over night

注意事项 (Note)

For research use only .