

产品名称: RAB43 Rabbit Polyclonal Antibody
产品货号: AP Rab16790



产品概述 (Summary)

产品名称 (Production Name)	RAB43 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAB43 RAB41
别名 (Alternative Names)	
基因 ID (Gene ID)	339122.0
蛋白 ID (SwissProt ID)	Q86YS6.Synthesized peptide derived from human protein . at AA range: 90-170

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	23kDa

研究背景 (Background)

similarity:Belongs to the small GTPase superfamily. Rab family.,tissue specificity:Widely expressed in brain, testis, lung,

产品名称: RAB43 Rabbit Polyclonal Antibody
产品货号: AP Rab16790




heart, ovary, colon, kidney, uterus and spleen but not in liver.,similarity:Belongs to the small GTPase superfamily. Rab family.,tissue specificity:Widely expressed in brain, testis, lung, heart, ovary, colon, kidney, uterus and spleen but not in liver.,

研究领域 (Research Area)

Signaling Pathway; G Protein Signaling; Small G Proteins; Ras Family; Signal Transduction; Protein Trafficking; Nuclear Import / Export; Neuroscience; Processes

图片 (Image Data)



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4°over night

注意事项 (Note)

For research use only .