

产品名称: Rabaptin-5 β Rabbit Polyclonal Antibody
产品货号: APRab16810



产品概述 (Summary)

产品名称 (Production Name)	Rabaptin-5 β Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RABEP2
别名 (Alternative Names)	RABEP2; RABPT5B; Rab GTPase-binding effector protein 2; Rabaptin-5beta
基因 ID (Gene ID)	79874.0
蛋白 ID (SwissProt ID)	Q9H5N1.The antiserum was produced against synthesized peptide derived from human RABEP2. AA range:353-402

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	64kDa

研究背景 (Background)


产品名称: Rabaptin-5 β Rabbit Polyclonal Antibody
产品货号: AP Rab16810




function: Plays a role in membrane trafficking and in homotypic early endosome fusion., similarity: Belongs to the rabaptin family., subunit: Heterodimer with RABGEF1. The dimer binds RAB5A that has been activated by GTP-binding., function: Plays a role in membrane trafficking and in homotypic early endosome fusion., similarity: Belongs to the rabaptin family., subunit: Heterodimer with RABGEF1. The dimer binds RAB5A that has been activated by GTP-binding.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from HeLa, Jurkat, 293, and HUVEC cells, using RABEP2 Antibody. The lane on the right is blocked with the synthesized peptide.




Western blot analysis of the lysates from HeLa cells using RABEP2 antibody.

产品名称: Rabaptin-5 β Rabbit Polyclonal Antibody
产品货号: APRab16810



Western Blot analysis of various cells using Rabaptin-5 β Polyclonal Antibody



Immunohistochemical analysis of paraffin-embedded Human breast cancer. Antibody was diluted at 1:100 (4 $^{\circ}$,overnight) .
High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negative contrl (right) obtaned from
antibody was pre-absorbed by immunogen peptide.

注意事项 (Note)

For research use only .