

产品名称: Rabr Rabbit Polyclonal Antibody
产品货号: AP Rab16819



产品概述 (Summary)

产品名称 (Production Name)	Rabr Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AGFG2 AGFG2; HRBL; RABR; Arf-GAP domain and FG repeat-containing protein 2;
别名 (Alternative Names)	HIV-1 Rev-binding protein-like protein; Rev/Rex activation domain-binding protein related; RAB-R
基因 ID (Gene ID)	3268.0
蛋白 ID (SwissProt ID)	O95081.The antiserum was produced against synthesized peptide derived from human AGFG2. AA range:91-140

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	49kDa

产品名称: Rabr Rabbit Polyclonal Antibody
产品货号: AP Rab16819




研究背景 (Background)

ArfGAP with FG repeats 2 (AGFG2) Homo sapiens This gene is a member of the HIV-1 Rev binding protein (HRB) family and encodes a protein with one Arf-GAP zinc finger domain, several phe-gly (FG) motifs, and four asn-pro-phe (NPF) motifs. This protein interacts with Eps15 homology (EH) domains and plays a role in the Rev export pathway, which mediates the nucleocytoplasmic transfer of proteins and RNAs. Alternatively spliced transcript variants have been described, but their biological validity has not been determined. [provided by RefSeq, Feb 2013], domain: Contains FG repeats and 4 N-P-F repeats., similarity: Contains 1 Arf-GAP domain., subunit: Interacts with EPS15R.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from MCF-7 cells, using AGFG2 Antibody. The lane on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .