

产品名称: Raftlin Rabbit Polyclonal Antibody
产品货号: APRab16858



产品概述 (Summary)

产品名称 (Production Name)	Raftlin Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RFTN1
别名 (Alternative Names)	RFTN1; KIAA0084; MIG2; Raftlin; Cell migration-inducing gene 2 protein; Raft-linking protein
基因 ID (Gene ID)	23180.0
蛋白 ID (SwissProt ID)	Q14699.The antiserum was produced against synthesized peptide derived from human RFTN1. AA range:432-481

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	65kDa

产品名称: Raftlin Rabbit Polyclonal Antibody
产品货号: APRab16858




研究背景 (Background)

function: May play a pivotal role in the formation and/or maintenance of lipid rafts. May regulate B-cell antigen receptor-mediated signaling., similarity: Belongs to the raftlin family., subcellular location: Associates with lipid rafts., tissue specificity: Specifically expressed by B-cells (at protein level)., function: May play a pivotal role in the formation and/or maintenance of lipid rafts. May regulate B-cell antigen receptor-mediated signaling., similarity: Belongs to the raftlin family., subcellular location: Associates with lipid rafts., tissue specificity: Specifically expressed by B-cells (at protein level).,


研究领域 (Research Area)

Tags & Cell Markers; Subcellular Markers; Plasma Membrane

图片 (Image Data)




Western blot analysis of the lysates from RAW264.7 cells using RFTN1 antibody.



Western blot analysis of the lysates from RAW264.7 cells using RFTN1 antibody.

产品名称: Raftlin Rabbit Polyclonal Antibody
产品货号: APRab16858



Western Blot analysis of 3T3 cells using Raftlin Polyclonal Antibody

注意事项 (Note)

For research use only .