

产品名称: Rap 2C Rabbit Polyclonal Antibody
产品货号: APRab16890



产品概述 (Summary)

产品名称 (Production Name)	Rap 2C Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAP2C
别名 (Alternative Names)	RAP2C; Ras-related protein Rap-2c
基因 ID (Gene ID)	57826.0
蛋白 ID (SwissProt ID)	Q9Y3L5.The antiserum was produced against synthesized peptide derived from human RAP2C. AA range:97-146

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	26kDa

研究背景 (Background)

产品名称: Rap 2C Rabbit Polyclonal Antibody
产品货号: AP Rab16890





The protein encoded by this gene is a member of the Ras-related protein subfamily of the Ras GTPase superfamily. Members of this family are small GTPases that act as molecular switches to regulate cellular proliferation, differentiation, and apoptosis. This protein has been reported to activate in vitro transcriptional activity of the serum response element. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2012], similarity: Belongs to the small GTPase superfamily. Ras family., tissue specificity: Expressed in mononuclear leukocytes and megakaryocytes. Not detected in platelets.,


研究领域 (Research Area)

Signal Transduction; Signaling Pathway; G Protein Signaling; Small G Proteins; Ras Family

图片 (Image Data)



产品名称: Rap 2C Rabbit Polyclonal Antibody
产品货号: AP Rab16890



Western Blot analysis of various cells using Rap 2C Polyclonal Antibody diluted at 1: 500

注意事项 (Note)

For research use only .