

产品名称: RASL10B Rabbit Polyclonal Antibody
产品货号: APRab16916



产品概述 (Summary)

产品名称 (Production Name)	RASL10B Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RASL10B
别名 (Alternative Names)	RASL10B; Ras-like protein family member 10B; Ras-like protein VTS58635; Ras-related protein 17; RRP17
基因 ID (Gene ID)	91608.0
蛋白 ID (SwissProt ID)	Q96S79.Synthesized peptide derived from RASL10B . at AA range: 80-160

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	22kDa

研究背景 (Background)


产品名称: RASL10B Rabbit Polyclonal Antibody
产品货号: APRab16916



similarity:Belongs to the small GTPase superfamily. Ras family.,similarity:Belongs to the small GTPase superfamily. Ras family.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of RASL10B Antibody. The lane on the right is blocked with the RASL10B peptide.

注意事项 (Note)

For research use only .