

产品名称: RASSF4 Rabbit Polyclonal Antibody
产品货号: APRab16923



产品概述 (Summary)

产品名称 (Production Name)	RASSF4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RASSF4
别名 (Alternative Names)	RASSF4; AD037; Ras association domain-containing protein 4
基因 ID (Gene ID)	83937.0
蛋白 ID (SwissProt ID)	Q9H2L5.The antiserum was produced against synthesized peptide derived from human RASSF4. AA range:134-183

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	37kDa

研究背景 (Background)


产品名称: RASSF4 Rabbit Polyclonal Antibody
产品货号: AP Rab16923




The function of this gene has not yet been determined but may involve a role in tumor suppression. Alternative splicing of this gene results in several transcript variants; however, most of the variants have not been fully described. [provided by RefSeq, Jul 2008],function:Potential tumor suppressor. May act as a KRAS effector protein. May promote apoptosis and cell cycle arrest.,similarity:Contains 1 Ras-associating domain.,similarity:Contains 1 SARAH domain.,subunit:Interacts directly with activated KRAS in a GTP-dependent manner.,tissue specificity:Widely expressed. Frequently down-regulated in tumor cell lines.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from COLO and HUVEC cells, using RASSF4 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from RAW264.7 cells using RASSF4 antibody.

产品名称: RASSF4 Rabbit Polyclonal Antibody
产品货号: APRab16923



Western Blot analysis of various cells using RASSF4 Polyclonal Antibody

注意事项 (Note)

For research use only .