

产品名称: RASSF6 Rabbit Polyclonal Antibody
产品货号: APRab16924



产品概述 (Summary)

产品名称 (Production Name)	RASSF6 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RASSF6
别名 (Alternative Names)	RASSF6; Ras association domain-containing protein 6
基因 ID (Gene ID)	166824.0
蛋白 ID (SwissProt ID)	Q6ZTQ3.The antiserum was produced against synthesized peptide derived from human RASSF6. AA range:129-178

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	45kDa

研究背景 (Background)


产品名称: RASSF6 Rabbit Polyclonal Antibody
产品货号: APRab16924




This gene encodes a member of the Ras-association domain family (RASSF). Members of this family form the core of a highly conserved tumor suppressor network, the Salvador-Warts-Hippo (SWH) pathway. The protein encoded by this gene is a Ras effector protein that induces apoptosis. A genomic region containing this gene has been linked to susceptibility to viral bronchiolitis. Alternative splicing results in multiple transcript variants and protein isoforms. [provided by RefSeq, Jul 2012],function:May act as a Ras effector protein.,sequence caution:Translated as Gln.,similarity:Contains 1 Ras-associating domain.,similarity:Contains 1 SARAH domain.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from HeLa cells, using RASSF6 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from K562 cells using RASSF6 antibody.

产品名称: RASSF6 Rabbit Polyclonal Antibody
产品货号: APRab16924



Western Blot analysis of various cells using RASSF6 Polyclonal Antibody

注意事项 (Note)

For research use only .