

产品名称: RASSF7 Rabbit Polyclonal Antibody
产品货号: APRab16925



产品概述 (Summary)

产品名称 (Production Name)	RASSF7 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RASSF7
别名 (Alternative Names)	RASSF7; C11orf13; HRC1; Ras association domain-containing protein 7; HRAS1-related cluster protein 1
基因 ID (Gene ID)	8045.0
蛋白 ID (SwissProt ID)	Q02833.The antiserum was produced against synthesized peptide derived from human RASSF7. AA range:159-208

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	

产品名称: RASSF7 Rabbit Polyclonal Antibody
产品货号: AP Rab16925




研究背景 (Background)

similarity:Contains 1 Ras-associating domain.,similarity:Contains 1 Ras-associating domain.,

研究领域 (Research Area)

Cell Biology; Cell Cycle; Cell Division; Spindle; Tags & Cell Markers; Subcellular Markers; Organelles; Centrosome; Epigenetics and Nuclear Signaling; Chromosome Structure; Chromatin assembly

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using RASSF7 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .