

**产品名称: RBM16 Rabbit Polyclonal Antibody**  
**产品货号: APRab16950**



## 产品概述 (Summary)

|                        |                                  |
|------------------------|----------------------------------|
| 产品名称 (Production Name) | RBM16 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody       |
| 宿主 (Host)              | Rabbit                           |
| 应用 (Application)       | WB,IHC                           |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                  |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                                                                 |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | SCAF8<br>SCAF8; CCAP7; KIAA1116; RBM16; Protein SCAF8; CDC5L complex-associated protein 7; RNA-binding motif protein 16; SR-related and CTD-associated factor 8 |
| 别名 (Alternative Names) |                                                                                                                                                                 |
| 基因 ID (Gene ID)        | 22828.0                                                                                                                                                         |
| 蛋白 ID (SwissProt ID)   | Q9UPN6.Synthesized peptide derived from RBM16 . at AA range: 550-630                                                                                            |

## 产品应用 (Application)

|                          |                                |
|--------------------------|--------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,IHC 1:50-1:300 |
| 蛋白分子量 (Molecular Weight) | 140kDa                         |

产品名称: RBM16 Rabbit Polyclonal Antibody  
产品货号: APRab16950




## 研究背景 (Background)

PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 CID domain.,similarity:Contains 1 RRM (RNA recognition motif) domain.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 CID domain.,similarity:Contains 1 RRM (RNA recognition motif) domain.,

## 研究领域 (Research Area)

## 图片 (Image Data)



Western Blot analysis of JK cells using RBM16 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA) .

## 注意事项 (Note)

For research use only .