

产品名称: RCL Rabbit Polyclonal Antibody
产品货号: APRab16982



产品概述 (Summary)

产品名称 (Production Name)	RCL Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RCL
别名 (Alternative Names)	RCL; C6orf108; Deoxyribonucleoside 5'-monophosphate N-glycosidase; c-Myc-responsive protein Rcl
基因 ID (Gene ID)	10591.0
蛋白 ID (SwissProt ID)	O43598.The antiserum was produced against synthesized peptide derived from human RCL. AA range:1-50

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	22kDa

产品名称: RCL Rabbit Polyclonal Antibody
产品货号: APRab16982




研究背景 (Background)

This gene was identified on the basis of its stimulation by c-Myc protein. The latter is a transcription factor that participates in the regulation of cell proliferation, differentiation, and apoptosis. The exact function of this gene is not known but studies in rat suggest a role in cellular proliferation and c-Myc-mediated transformation. Two alternative transcripts encoding different proteins have been described. [provided by RefSeq, Jul 2008],caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,function:May play a role during cellular proliferation and c-Myc-mediated transformation.,tissue specificity:Expression is elevated in lymphoid cells that overexpress c-Myc.,


研究领域 (Research Area)

Other Antibodies; Cell Biology; Cell Cycle; Cell Division; Other cell division antibodies

图片 (Image Data)




Western blot analysis of lysates from Jurkat and K562 cells, using RCL Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from 293 cells using RCL antibody.

产品名称: RCL Rabbit Polyclonal Antibody
产品货号: APRab16982



Western Blot analysis of various cells using RCL Polyclonal Antibody

注意事项 (Note)

For research use only .