

产品名称: RECS1 Rabbit Polyclonal Antibody
产品货号: APRab16994



产品概述 (Summary)

产品名称 (Production Name)	RECS1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TMBIM1
别名 (Alternative Names)	TMBIM1; LFG3; RECS1; PP1201; PSEC0158; Protein lifeguard 3; Protein RECS1 homolog; Transmembrane BAX inhibitor motif-containing protein 1
基因 ID (Gene ID)	64114.0
蛋白 ID (SwissProt ID)	Q969X1.The antiserum was produced against synthesized peptide derived from human TMBIM1. AA range:251-300

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	35kDa

产品名称: RECS1 Rabbit Polyclonal Antibody
产品货号: APRab16994




研究背景 (Background)

similarity:Belongs to the BI1 family.,similarity:Belongs to the BI1 family.,


研究领域 (Research Area)

Cell Biology; Other Antibodies


图片 (Image Data)



Western blot analysis of lysates from HUVEC cells, using TMBIM1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from K562 cells using TMBIM1 antibody.



Western Blot analysis of various cells using RECS1 Polyclonal Antibody diluted at 1: 1000

产品名称: RECS1 Rabbit Polyclonal Antibody
产品货号: APRab16994



注意事项 (Note)

For research use only .