

产品名称: REXO1 Rabbit Polyclonal Antibody
产品货号: APRab17042



产品概述 (Summary)

产品名称 (Production Name)	REXO1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	REXO1 REXO1; ELOABP1; KIAA1138; TCEB3BP1; RNA exonuclease 1 homolog;
别名 (Alternative Names)	Elongin-A-binding protein 1; EloA-BP1; Transcription elongation factor B polypeptide 3-binding protein 1
基因 ID (Gene ID)	57455.0
蛋白 ID (SwissProt ID)	Q8N1G1. The antiserum was produced against synthesized peptide derived from human REXO1. AA range: 185-234

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

产品名称: REXO1 Rabbit Polyclonal Antibody
产品货号: APRab17042




研究背景 (Background)

function: Seems to have no detectable effect on transcription elongation in vitro., similarity: Belongs to the REXO1/REXO3 family., similarity: Contains 1 exonuclease domain., subunit: Interacts with TCEA2 and TCEB3., tissue specificity: Ubiquitously expressed., function: Seems to have no detectable effect on transcription elongation in vitro., similarity: Belongs to the REXO1/REXO3 family., similarity: Contains 1 exonuclease domain., subunit: Interacts with TCEA2 and TCEB3., tissue specificity: Ubiquitously expressed.,

研究领域 (Research Area)

Epigenetics and Nuclear Signaling; Transcription; Elongation factors

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human brain, using REXO1 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .