

产品名称: RFPL4A Rabbit Polyclonal Antibody
产品货号: APRab17059



产品概述 (Summary)

产品名称 (Production Name)	RFPL4A Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RFPL4A
别名 (Alternative Names)	RFPL4A; RFPL4; RNF210; Ret finger protein-like 4A; RING finger protein 210
基因 ID (Gene ID)	342931.0
蛋白 ID (SwissProt ID)	A6NLU0.The antiserum was produced against synthesized peptide derived from human RFPL4A. AA range:155-204

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	42kDa

研究背景 (Background)


产品名称: RFPL4A Rabbit Polyclonal Antibody
产品货号: APRab17059




The ret finger protein-like 4 gene, Rfpl4, encodes a putative E3 ubiquitin-protein ligase expressed in adult germ cells.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of lysates from HUVEC cells, using RFPL4A Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HeLa cells using RFPL4A antibody.



Western Blot analysis of various cells using RFPL4A Polyclonal Antibody

产品名称: RFPL4A Rabbit Polyclonal Antibody
产品货号: APRab17059



注意事项 (Note)

For research use only .