

产品名称: RGAG1 Rabbit Polyclonal Antibody
产品货号: APRab17070



产品概述 (Summary)

产品名称 (Production Name)	RGAG1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RGAG1
别名 (Alternative Names)	RGAG1; KIAA1318; Retrotransposon gag domain-containing protein 1; Tumor antigen BJ-HCC-23
基因 ID (Gene ID)	57529.0
蛋白 ID (SwissProt ID)	Q8NET4.The antiserum was produced against synthesized peptide derived from human RGAG1. AA range:1331-1380

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	144kDa

产品名称: RGAG1 Rabbit Polyclonal Antibody
产品货号: APRab17070




研究背景 (Background)


RGAG1 (Retrotransposon Gag Domain Containing 1) is a Protein Coding gene. Diseases associated with RGAG1 include non-syndromic x-linked intellectual disability.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of lysates from COLO cells, using RGAG1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HT-29 cells using RGAG1 antibody.



产品名称: RGAG1 Rabbit Polyclonal Antibody
产品货号: APRab17070



Western Blot analysis of various cells using RGAG1 Polyclonal Antibody

注意事项 (Note)

For research use only .