

产品名称: RGP2/3/4/6/8 Rabbit Polyclonal Antibody
产品货号: APRab17084



产品概述 (Summary)

产品名称 (Production Name)	RGP2/3/4/6/8 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RGP2
别名 (Alternative Names)	
基因 ID (Gene ID)	729540/84220/729857 O14715/Q53T03/Q7Z3J3/A6NKT7/P0C839.The antiserum was produced against synthesized peptide derived from human RGP2/3/4/6/8. AA range:577-626
蛋白 ID (SwissProt ID)	

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	110kDa

产品名称: RGP2/3/4/6/8 Rabbit Polyclonal Antibody
产品货号: APRab17084




研究背景 (Background)


caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,miscellaneous:One of the 8 copies of RANBP2 clustered close to the centromere of chromosome 2.,miscellaneous:One of the 8 copies of RANBP2 clustered close to the chromosome 2 centromere.,similarity:Contains 1 GRIP domain.,similarity:Contains 2 RanBD1 domains.,similarity:Contains 2 TPR repeats.,similarity:Contains 3 TPR repeats.,tissue specificity:Expressed in testis.,caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,miscellaneous:One of the 8 copies of RANBP2 clustered close to the centromere of chromosome 2.,miscellaneous:One of the 8 copies of RANBP2 clustered close to the chromosome 2 centromere.,similarity:Contains 1 GRIP domain.,similarity:Contains 2 RanBD1 domains.,similarity:Contains 2 TPR repeats.,similarity:Contains 3 TPR repeats.,tissue specificity:Expressed in testis.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of RGP2/3/4/6/8 Antibody. The lane on the right is blocked with the RGP2/3/4/6/8 peptide.



Western blot analysis of the lysates from Jurkat cells using RGP2/3/4/6/8 antibody.

产品名称: RGP2/3/4/6/8 Rabbit Polyclonal Antibody
产品货号: APRab17084



注意事项 (Note)

For research use only .