

产品名称: Rho H Rabbit Polyclonal Antibody
产品货号: APRab17126



产品概述 (Summary)

产品名称 (Production Name)	Rho H Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RhoH
别名 (Alternative Names)	RHOH; ARHH; TTF; Rho-related GTP-binding protein RhoH; GTP-binding protein TTF; Translocation three four protein
基因 ID (Gene ID)	399.0
蛋白 ID (SwissProt ID)	Q15669.The antiserum was produced against synthesized peptide derived from human RhoH. AA range:141-190

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	21kDa

产品名称: Rho H Rabbit Polyclonal Antibody
产品货号: APRab17126




研究背景 (Background)

The protein encoded by this gene is a member of the Ras superfamily of guanosine triphosphate (GTP)-metabolizing enzymes. The encoded protein is expressed in hematopoietic cells, where it functions as a negative regulator of cell growth and survival. This gene may be hypermutated or misexpressed in leukemias and lymphomas. Chromosomal translocations in non-Hodgkin's lymphoma occur between this locus and B-cell CLL/lymphoma 6 (BCL6) on chromosome 3, leading to the production of fusion transcripts. Alternative splicing in the 5' untranslated region results in multiple transcript variants that encode the same protein. [provided by RefSeq, May 2013], disease: A chromosomal aberration involving RHOH is found in a non-Hodgkin lymphoma cell line. Translocation t(3;4)(q27;p11) with BCL6, similarity: Belongs to the small GTPase superfamily. Rho family, tissue specificity: Transcribed only in hemopoietic cells,


研究领域 (Research Area)

Leukocyte transendothelial migration;

图片 (Image Data)




Immunofluorescence analysis of A549 cells, using RhoH Antibody. The picture on the right is blocked with the synthesized



Western blot analysis of lysates from HT-29 cells, using RhoH Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: Rho H Rabbit Polyclonal Antibody
产品货号: APRab17126



Western blot analysis of the lysates from RAW264.7 cells using RhoH antibody.

注意事项 (Note)

For research use only .