

产品名称: RHOBTB3 Rabbit Polyclonal Antibody
产品货号: APRab17128



产品概述 (Summary)

产品名称 (Production Name)	RHOBTB3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RHOBTB3
别名 (Alternative Names)	RHOBTB3; KIAA0878; Rho-related BTB domain-containing protein 3
基因 ID (Gene ID)	22836.0
蛋白 ID (SwissProt ID)	O94955.The antiserum was produced against synthesized peptide derived from human RHOBTB3. AA range:335-384

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	70kDa

研究背景 (Background)


产品名称: RHOBTB3 Rabbit Polyclonal Antibody
产品货号: APRab17128




RHOBTB3 is a member of the evolutionarily conserved RHOBTB subfamily of Rho GTPases. For background information on RHOBTBs, see RHOBTB1 (MIM 607351).[supplied by OMIM, Apr 2004],similarity:Contains 2 BTB (POZ) domains.,tissue specificity:Ubiquitous. Highly expressed in neural and cardiac tissues, pancreas, placenta and testis.,

研究领域 (Research Area)


图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using RHOBTB3 Antibody. The picture on the right is blocked with the synthesized peptide.



Western Blot analysis of K562 cells using RHOBTB3 Polyclonal Antibody diluted at 1: 2000



产品名称: RHOBTB3 Rabbit Polyclonal Antibody
产品货号: APRab17128



Western blot analysis of customer 's lysis using RHOBTB3 antibody. Antibody was diluted at 1:2000

注意事项 (Note)

For research use only .