

产品名称: Ribosomal Protein L10L Rabbit Polyclonal Antibody
产品货号: APRab17144



产品概述 (Summary)

产品名称 (Production Name)	Ribosomal Protein L10L Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RPL10L
别名 (Alternative Names)	RPL10L; 60S ribosomal protein L10-like
基因 ID (Gene ID)	140801.0
蛋白 ID (SwissProt ID)	Q96L21.The antiserum was produced against synthesized peptide derived from human RPL10L. AA range:151-200

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	20kDa

产品名称: Ribosomal Protein L10L Rabbit Polyclonal Antibody
产品货号: APRab17144




研究背景 (Background)

This gene encodes a protein sharing sequence similarity with ribosomal protein L10. It is not currently known whether the encoded protein is a functional ribosomal protein or whether it has evolved a function that is independent of the ribosome. This gene is intronless. [provided by RefSeq, Jul 2008],function:May play a role in compensating for the inactivated X-linked gene during spermatogenesis.,miscellaneous:This gene has no introns in its coding regions, and therefore was most likely produced by retrotransposition of the original X-linked gene during evolution.,similarity:Belongs to the ribosomal protein L10e family.,tissue specificity:Testis specific.,


研究领域 (Research Area)

Ribosome;

图片 (Image Data)



Western blot analysis of lysates from Jurkat, COS7, and HeLa cells, using RPL10L Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HUVEC cells using RPL10L antibody.

产品名称: Ribosomal Protein L10L Rabbit Polyclonal
Antibody
产品货号: APRab17144



Western Blot analysis of various cells using Ribosomal Protein L10L Polyclonal Antibody diluted at 1 : 1000

注意事项 (Note)

For research use only .