

产品名称: Ribosomal Protein L26L Rabbit Polyclonal Antibody
产品货号: APRab17153



产品概述 (Summary)

产品名称 (Production Name)	Ribosomal Protein L26L Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RPL26L1
别名 (Alternative Names)	RPL26L1; RPL26P1; 60S ribosomal protein L26-like 1
基因 ID (Gene ID)	51121.0
蛋白 ID (SwissProt ID)	Q9UNX3.The antiserum was produced against synthesized peptide derived from human RPL26L. AA range:96-145

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	20kDa

产品名称: Ribosomal Protein L26L Rabbit Polyclonal Antibody
产品货号: APRab17153




研究背景 (Background)

This gene encodes a protein that shares high sequence similarity with ribosomal protein L26. Alternative splicing results in multiple transcript variants encoding the same protein. [provided by RefSeq, Dec 2015],similarity:Belongs to the ribosomal protein L24P family.,

研究领域 (Research Area)

Ribosome;

图片 (Image Data)




Western blot analysis of lysates from HUVEC and MCF-7 cells, using RPL26L Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using Ribosomal Protein L26L Polyclonal Antibody diluted at 1 : 2000

产品名称: Ribosomal Protein L26L Rabbit Polyclonal Antibody
产品货号: APRab17153



Immunohistochemical analysis of paraffin-embedded Human brain. Antibody was diluted at 1:100 (4°,overnight) . High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negative contrl (right) obtaned from antibody was pre-absorbed by immunogen peptide.

注意事项 (Note)

For research use only .