

产品名称: Ribosomal Protein L39L Rabbit Polyclonal Antibody
产品货号: APRab17165



产品概述 (Summary)

产品名称 (Production Name)	Ribosomal Protein L39L Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RPL39L
别名 (Alternative Names)	RPL39L; RPL39L1; 60S ribosomal protein L39-like; 60S ribosomal protein L39-2
基因 ID (Gene ID)	116832.0
蛋白 ID (SwissProt ID)	Q96EH5. The antiserum was produced against synthesized peptide derived from human RPL39L. AA range: 1-50

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	

产品名称: Ribosomal Protein L39L Rabbit Polyclonal Antibody
产品货号: APRab17165




研究背景 (Background)

This gene encodes a protein sharing high sequence similarity with ribosomal protein L39. Although the name of this gene has been referred to as 'ribosomal protein L39' in the public databases, its official name is 'ribosomal protein L39-like'. It is not currently known whether the encoded protein is a functional ribosomal protein or whether it has evolved a function that is independent of the ribosome. [provided by RefSeq, Jul 2008],similarity:Belongs to the ribosomal protein L39e family.,tissue specificity:Testis specific.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human brain, using RPL39L Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .