

产品名称: Ribosomal Protein L3L Rabbit Polyclonal Antibody
产品货号: APRab17166



产品概述 (Summary)

产品名称 (Production Name)	Ribosomal Protein L3L Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RPL3L
别名 (Alternative Names)	RPL3L; 60S ribosomal protein L3-like
基因 ID (Gene ID)	6123.0
蛋白 ID (SwissProt ID)	Q92901. The antiserum was produced against synthesized peptide derived from human RPL3L. AA range: 211-260

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:200-1:1000, ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)


产品名称: Ribosomal Protein L3L Rabbit Polyclonal Antibody
产品货号: APRab17166

This gene encodes a protein that shares sequence similarity with ribosomal protein L3. The protein belongs to the L3P family of ribosomal proteins. Unlike the ubiquitous expression of ribosomal protein genes, this gene has a tissue-specific pattern of expression, with the highest levels of expression in skeletal muscle and heart. It is not currently known whether the encoded protein is a functional ribosomal protein or whether it has evolved a function that is independent of the ribosome. [provided by RefSeq, Jul 2008],similarity:Belongs to the ribosomal protein L3P family.,

研究领域 (Research Area)

Ribosome;

图片 (Image Data)



Immunofluorescence analysis of MCF7 cells, using RPL3L Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .