

产品名称: Ribosomal Protein S19BP Rabbit Polyclonal Antibody
产品货号: APRab17178



产品概述 (Summary)

产品名称 (Production Name)	Ribosomal Protein S19BP Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RPS19BP1
别名 (Alternative Names)	RPS19BP1; AROS; Active regulator of SIRT1; 40S ribosomal protein S19-binding protein 1; RPS19-binding protein 1; S19BP
基因 ID (Gene ID)	91582.0
蛋白 ID (SwissProt ID)	Q86WX3. The antiserum was produced against synthesized peptide derived from human RPS19BP1. AA range: 71-120

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	

产品名称: Ribosomal Protein S19BP Rabbit Polyclonal Antibody
产品货号: AP Rab17178




研究背景 (Background)

function: Direct regulator of SIRT1. Enhances SIRT1-mediated deacetylation of p53/TP53, thereby participating in inhibition of p53/TP53-mediated transcriptional activity., similarity: Belongs to the AROS family., subunit: Interacts with RPS19 (By similarity). Interacts with SIRT1., tissue specificity: Widely expressed (at protein level)., function: Direct regulator of SIRT1. Enhances SIRT1-mediated deacetylation of p53/TP53, thereby participating in inhibition of p53/TP53-mediated transcriptional activity., similarity: Belongs to the AROS family., subunit: Interacts with RPS19 (By similarity). Interacts with SIRT1., tissue specificity: Widely expressed (at protein level).,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using RPS19BP1 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .