

**产品名称: Ribosomal Protein S27L Rabbit Polyclonal Antibody**  
**产品货号: APRab17184**



## 产品概述 (Summary)

|                        |                                                   |
|------------------------|---------------------------------------------------|
| 产品名称 (Production Name) | Ribosomal Protein S27L Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody                        |
| 宿主 (Host)              | Rabbit                                            |
| 应用 (Application)       | WB,IHC,ICC/IF,ELISA                               |
| 种属反应性 (Reactivity)     | Human,Mouse                                       |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                        |
|------------------------|--------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | RPS27L                                                                                                 |
| 别名 (Alternative Names) | RPS27L; 40S ribosomal protein S27-like                                                                 |
| 基因 ID (Gene ID)        | 51065.0                                                                                                |
| 蛋白 ID (SwissProt ID)   | Q71UM5.The antiserum was produced against synthesized peptide derived from human RPS27L. AA range:1-50 |

## 产品应用 (Application)

|                          |                                                                        |
|--------------------------|------------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000 |
| 蛋白分子量 (Molecular Weight) | 18kDa                                                                  |

产品名称: Ribosomal Protein S27L Rabbit Polyclonal Antibody  
产品货号: APRab17184




## 研究背景 (Background)

This gene encodes a protein sharing 96% amino acid similarity with ribosomal protein S27, which suggests the encoded protein may be a component of the 40S ribosomal subunit. [provided by RefSeq, Jul 2008],cofactor: Binds 1 zinc ion per subunit.,similarity: Belongs to the ribosomal protein S27e family.,


## 研究领域 (Research Area)

Ribosome;

## 图片 (Image Data)




Western blot analysis of lysates from COLO cells, using RPS27L Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from COLO205 cells using RPS27L antibody.

产品名称: Ribosomal Protein S27L Rabbit Polyclonal  
Antibody  
产品货号: AP Rab17184



Western Blot analysis of various cells using Ribosomal Protein S27L Polyclonal Antibody

#### 注意事项 (Note)

For research use only .