

产品名称: Ribosomal Protein S4X Rabbit Polyclonal Antibody
产品货号: APRab17186



产品概述 (Summary)

产品名称 (Production Name)	Ribosomal Protein S4X Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat,Cat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RPS4X
别名 (Alternative Names)	RPS4X; CCG2; RPS4; SCAR; 40S ribosomal protein S4; X isoform; SCR10; Single copy abundant mRNA protein
基因 ID (Gene ID)	6191.0
蛋白 ID (SwissProt ID)	P62701.The antiserum was produced against synthesized peptide derived from human RPS4X. AA range:81-130

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	30kDa

产品名称: Ribosomal Protein S4X Rabbit Polyclonal Antibody
产品货号: APRab17186




研究背景 (Background)

Cytoplasmic ribosomes, organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes ribosomal protein S4, a component of the 40S subunit. Ribosomal protein S4 is the only ribosomal protein known to be encoded by more than one gene, namely this gene and ribosomal protein S4, Y-linked (RPS4Y). The 2 isoforms encoded by these genes are not identical, but are functionally equivalent. Ribosomal protein S4 belongs to the S4E family of ribosomal proteins. This gene is not subject to X-inactivation. It has been suggested that haploinsufficiency of the ribosomal protein S4 genes plays a role in Turner syndrome; however, this hypothesis is controversial. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudomiscellaneous: This isoform is encoded on the X chromosome., similarity: Belongs to the ribosomal protein S4e family., similarity: Contains 1 S4 RNA-binding domain.,


研究领域 (Research Area)

Ribosome;

图片 (Image Data)




Immunofluorescence analysis of HUVEC cells, using RPS4X Antibody. The picture on the right is blocked with the synthesized




产品名称: Ribosomal Protein S4X Rabbit Polyclonal
Antibody
产品货号: APRab17186




Western blot analysis of lysates from HeLa and COLO cells, using RPS4X Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HT-29 cells using RPS4X antibody.



Western Blot analysis of various cells using Ribosomal Protein S4X Polyclonal Antibody diluted at 1: 1000



Western Blot analysis of customer 's (cat sample) using Ribosomal Protein S4X Polyclonal Antibody diluted at 1: 1000

注意事项 (Note)

For research use only .