

产品名称: Rim4 Rabbit Polyclonal Antibody
产品货号: APRab17195



产品概述 (Summary)

产品名称 (Production Name)	Rim4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RIMS4
别名 (Alternative Names)	RIMS4; C20orf190; Regulating synaptic membrane exocytosis protein 4; RIM4 gamma; Rab3-interacting molecule 4; RIM 4
基因 ID (Gene ID)	140730.0
蛋白 ID (SwissProt ID)	Q9H426.The antiserum was produced against synthesized peptide derived from human RIMS4. AA range:33-82

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	30kDa

产品名称: Rim4 Rabbit Polyclonal Antibody
产品货号: APRab17195




研究背景 (Background)

caution:Does not bind Rab-3.,function:Regulates synaptic membrane exocytosis.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 C2 domain.,subunit:Binds PPPIA3.,caution:Does not bind Rab-3.,function:Regulates synaptic membrane exocytosis.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 C2 domain.,subunit:Binds PPPIA3.,

研究领域 (Research Area)

图片 (Image Data)




Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using RIMS4 Antibody. The picture on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using Rim4 Polyclonal Antibody diluted at 1 : 500

产品名称: Rim4 Rabbit Polyclonal Antibody
产品货号: AP Rab17195



Western blot analysis of 3T3 lysis using Rim4 antibody. Antibody was diluted at 1:500

注意事项 (Note)

For research use only .