

产品名称: RNF113B Rabbit Polyclonal Antibody
产品货号: APRab17284



产品概述 (Summary)

产品名称 (Production Name)	RNF113B Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RNF113B
别名 (Alternative Names)	RNF113B; RNF161; ZNF183L1; RING finger protein 113B; Zinc finger protein 183-like 1
基因 ID (Gene ID)	140432.0
蛋白 ID (SwissProt ID)	Q8IZP6.Synthesized peptide derived from the Internal region of human RNF113B.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	36kDa

产品名称: RNF113B Rabbit Polyclonal Antibody
产品货号: AP Rab17284




研究背景 (Background)


similarity:Contains 1 C3H1-type zinc finger.,similarity:Contains 1 RING-type zinc finger.,similarity:Contains 1 C3H1-type zinc finger.,similarity:Contains 1 RING-type zinc finger.,

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of various cells using RNF113B Polyclonal Antibody



Western Blot analysis of 293 cells using RNF113B Polyclonal Antibody

注意事项 (Note)

For research use only .