

产品名称: RNF130 Rabbit Polyclonal Antibody
产品货号: APRab17286



产品概述 (Summary)

产品名称 (Production Name)	RNF130 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RNF130
别名 (Alternative Names)	RNF130; E3 ubiquitin-protein ligase RNF130; Goliath homolog; H-Goliath; RING finger protein 130
基因 ID (Gene ID)	55819.0
蛋白 ID (SwissProt ID)	Q86XS8.The antiserum was produced against synthesized peptide derived from human RNF130. AA range:231-280

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	36kDa

产品名称: RNF130 Rabbit Polyclonal Antibody
产品货号: APRab17286




研究背景 (Background)


ring finger protein 130(RNF130) Homo sapiens The protein encoded by this gene contains a RING finger motif and is similar to g1, a Drosophila zinc-finger protein that is expressed in mesoderm and involved in embryonic development. The expression of the mouse counterpart was found to be upregulated in myeloblastic cells following IL3 deprivation, suggesting that this gene may regulate growth factor withdrawal-induced apoptosis of myeloid precursor cells. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013],function:May have a role during the programmed cell death of hematopoietic cells.,similarity:Contains 1 PA (protease associated) domain.,similarity:Contains 1 RING-type zinc finger.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from HeLa, Jurkat, HUVEC, and MCF-7 cells, using RNF130 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using RNF130 Polyclonal Antibody diluted at 1: 1000

注意事项 (Note)

For research use only .

产品名称: RNF130 Rabbit Polyclonal Antibody
产品货号: APRab17286

