

产品名称: RNF149 Rabbit Polyclonal Antibody
产品货号: APRab17289



产品概述 (Summary)

产品名称 (Production Name)	RNF149 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RNF149
别名 (Alternative Names)	RNF149; DNAPTP2; E3 ubiquitin-protein ligase RNF149; DNA polymerase-transactivated protein 2; RING finger protein 149
基因 ID (Gene ID)	284996.0
蛋白 ID (SwissProt ID)	Q8NC42. The antiserum was produced against synthesized peptide derived from human RNF149. AA range: 211-260

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	

产品名称: RNF149 Rabbit Polyclonal Antibody
产品货号: APRab17289




研究背景 (Background)

domain:The RING-type zinc finger domain mediates binding to an E2 ubiquitin-conjugating enzyme.,pathway:Protein modification; protein ubiquitination.,similarity:Contains 1 PA (protease associated) domain.,similarity:Contains 1 RING-type zinc finger.,domain:The RING-type zinc finger domain mediates binding to an E2 ubiquitin-conjugating enzyme.,pathway:Protein modification; protein ubiquitination.,similarity:Contains 1 PA (protease associated) domain.,similarity:Contains 1 RING-type zinc finger.,

研究领域 (Research Area)

Cell Biology; Proteolysis / Ubiquitin; Proteasome / Ubiquitin; Ubiquitin E3 Enzymes; RING Finger E3 Ligase

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using RNF149 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .