

产品名称: RNF15 Rabbit Polyclonal Antibody
产品货号: APRab17290



产品概述 (Summary)

产品名称 (Production Name)	RNF15 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TRIM38
别名 (Alternative Names)	TRIM38; RNF15; RORET; Tripartite motif-containing protein 38; RING finger protein 15; Zinc finger protein RoRet
基因 ID (Gene ID)	10475.0
蛋白 ID (SwissProt ID)	O00635.The antiserum was produced against synthesized peptide derived from human TRIM38. AA range:161-210

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	58kDa

产品名称: RNF15 Rabbit Polyclonal Antibody
产品货号: AP Rab17290





研究背景 (Background)

This gene encodes a member of the tripartite motif (TRIM) family. The encoded protein contains a RING-type zinc finger, B box-type zinc finger and SPRY domain. The function of this protein has not been identified. A pseudogene of this gene is located on the long arm of chromosome 4. [provided by RefSeq, Jul 2012], similarity: Contains 1 B box-type zinc finger., similarity: Contains 1 B30.2/SPRY domain., similarity: Contains 1 RING-type zinc finger., tissue specificity: Ubiquitous.,

研究领域 (Research Area)

图片 (Image Data)



注意事项 (Note)

For research use only .