

产品名称: RNF41 Rabbit Polyclonal Antibody
产品货号: AP Rab17294



产品概述 (Summary)

产品名称 (Production Name)	RNF41 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RNF41 FLRF NRDP1 SBBIO3
别名 (Alternative Names)	
基因 ID (Gene ID)	10193.0
蛋白 ID (SwissProt ID)	Q9H4P4.Synthesized peptide derived from part region of human protein

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	34kDa

研究背景 (Background)

This gene encodes an E3 ubiquitin ligase. The encoded protein plays a role in type 1 cytokine receptor signaling by controlling the balance between JAK2-associated cytokine receptor degradation and ectodomain shedding. Alternative

产品名称: RNF41 Rabbit Polyclonal Antibody
产品货号: AP Rab17294




splicing results in multiple transcript variants. [provided by RefSeq, Jun 2011],function:Acts as E3 ubiquitin-protein ligase and regulates the degradation of target proteins, such as ERBB3.,pathway:Protein modification; protein ubiquitination.,PTM:Autoubiquitinated.,similarity:Contains 1 RING-type zinc finger.,similarity:Contains 1 SIAH-type zinc finger.,subunit:Interacts with USP8 and ERBB3.,tissue specificity:Detected in ovary, testis and prostate.,

研究领域 (Research Area)

Endocytosis;

图片 (Image Data)



Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4°over night

注意事项 (Note)

For research use only .