

产品名称: RPA14 Rabbit Polyclonal Antibody
产品货号: AP Rab17336



产品概述 (Summary)

产品名称 (Production Name)	RPA14 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RPA3
别名 (Alternative Names)	Replication protein A 14 kDa subunit (RP-A p14) (Replication factor A protein 3) (RF-A protein 3)
基因 ID (Gene ID)	6119.0
蛋白 ID (SwissProt ID)	P35244.Synthetic peptide from human protein at AA range: 30-80

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	14kDa

研究背景 (Background)

产品名称: RPA14 Rabbit Polyclonal Antibody
产品货号: APRab17336




function:Required for DNA recombination, repair and replication. The activity of RP-A is mediated by single-stranded DNA binding and protein interactions.,subunit:Heterotrimer of 70, 32 and 14 kDa chains. The DNA-binding activity may reside exclusively on the 70 kDa subunit. Interacts with RPA4.,function:Required for DNA recombination, repair and replication. The activity of RP-A is mediated by single-stranded DNA binding and protein interactions.,subunit:Heterotrimer of 70, 32 and 14 kDa chains. The DNA-binding activity may reside exclusively on the 70 kDa subunit. Interacts with RPA4.,

研究领域 (Research Area)

DNA replication;Nucleotide excision repair;Mismatch repair;Homologous recombination;

图片 (Image Data)



Western blot analysis of lysate, antibody was diluted at 2000. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .