

产品名称: RPAP2 Rabbit Polyclonal Antibody
产品货号: APRab17346



产品概述 (Summary)

产品名称 (Production Name)	RPAP2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RPAP2
别名 (Alternative Names)	RPAP2; C1orf82; Putative RNA polymerase II subunit B1 CTD phosphatase RPAP2; RNA polymerase II-associated protein 2
基因 ID (Gene ID)	79871.0
蛋白 ID (SwissProt ID)	Q8IXW5.The antiserum was produced against synthesized peptide derived from human RPAP2. AA range:460-509

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	70kDa

产品名称: RPAP2 Rabbit Polyclonal Antibody
产品货号: APRab17346




研究背景 (Background)


function:Forms an interface between the RNA polymerase II enzyme and chaperone/scaffolding protein, suggesting that it is required to connect RNA polymerase II to regulators of protein complex formation.,similarity:Belongs to the RPAP2 family.,subunit:Tightly associated with the RNA polymerase II complex.,function:Forms an interface between the RNA polymerase II enzyme and chaperone/scaffolding protein, suggesting that it is required to connect RNA polymerase II to regulators of protein complex formation.,similarity:Belongs to the RPAP2 family.,subunit:Tightly associated with the RNA polymerase II complex.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of RPAP2 Antibody. The lane on the right is blocked with the RPAP2 peptide.



Western Blot analysis of various cells using RPAP2 Polyclonal Antibody

注意事项 (Note)

For research use only .