

产品名称: RRP22 Rabbit Polyclonal Antibody
产品货号: AP Rab17387



产品概述 (Summary)

产品名称 (Production Name)	RRP22 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RASL10A
别名 (Alternative Names)	RASL10A; RRP22; Ras-like protein family member 10A; Ras-like protein RRP22; Ras-related protein on chromosome 22
基因 ID (Gene ID)	10633.0
蛋白 ID (SwissProt ID)	Q92737.The antiserum was produced against synthesized peptide derived from human RASL10A. AA range:107-156

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	30kDa

产品名称: RRP22 Rabbit Polyclonal Antibody
产品货号: APRab17387




研究背景 (Background)


similarity:Belongs to the small GTPase superfamily. Ras family.,similarity:Belongs to the small GTPase superfamily. Ras family.,

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of lysates from HeLa cells, using RASL10A Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from 293 cells using RASL10A antibody.



产品名称: RRP22 Rabbit Polyclonal Antibody
产品货号: APRab17387



Western Blot analysis of MOUSE-SPLEEN cells using RRP22 Polyclonal Antibody diluted at 1 : 1000

注意事项 (Note)

For research use only .