

产品名称: RRS1 Rabbit Polyclonal Antibody
产品货号: APRab17391



产品概述 (Summary)

产品名称 (Production Name)	RRS1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RRS1
别名 (Alternative Names)	RRS1; KIAA0112; RRR; Ribosome biogenesis regulatory protein homolog
基因 ID (Gene ID)	23212.0
蛋白 ID (SwissProt ID)	Q15050.The antiserum was produced against synthesized peptide derived from human RRS1. AA range:254-303

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	41kDa

研究背景 (Background)

产品名称: RRS1 Rabbit Polyclonal Antibody
产品货号: APRab17391




function:Involved in ribosome biogenesis.,similarity:Belongs to the RRS1 family.,function:Involved in ribosome biogenesis.,similarity:Belongs to the RRS1 family.,


研究领域 (Research Area)

Ribosome; Epigenetics and Nuclear Signaling; DNA / RNA; Translation; Regulation


图片 (Image Data)



Western blot analysis of lysates from HUVEC and 293 cells, using RRS1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HT-29 cells using RRS1 antibody.




Western Blot analysis of various cells using RRS1 Polyclonal Antibody diluted at 1: 2000 cells nucleus extracted by Minute

产品名称: RRS1 Rabbit Polyclonal Antibody
产品货号: APRab17391



TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA) .



Western blot analysis of HELA KB SH-SY5Y lysis using RRS1 antibody. Antibody was diluted at 1:2000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA) .

注意事项 (Note)

For research use only .