

产品名称: RSAD1 Rabbit Polyclonal Antibody
产品货号: APRab17400



产品概述 (Summary)

产品名称 (Production Name)	RSAD1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RSAD1
别名 (Alternative Names)	RSAD1; Radical S-adenosyl methionine domain-containing protein 1; mitochondrial; Oxygen-independent coproporphyrinogen-III oxidase-like protein RSAD1
基因 ID (Gene ID)	55316.0
蛋白 ID (SwissProt ID)	Q9HA92.The antiserum was produced against synthesized peptide derived from human RSAD1. AA range:298-347

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	52kDa

产品名称: RSAD1 Rabbit Polyclonal Antibody
产品货号: APRab17400




研究背景 (Background)


cofactor: Binds 1 4Fe-4S cluster. The cluster is coordinated with 3 cysteines and an exchangeable S-adenosyl-L-methionine .,function: May be involved in porphyrin cofactor biosynthesis.,similarity: Belongs to the anaerobic coproporphyrinogen-III oxidase family.,cofactor: Binds 1 4Fe-4S cluster. The cluster is coordinated with 3 cysteines and an exchangeable S-adenosyl-L-methionine .,function: May be involved in porphyrin cofactor biosynthesis.,similarity: Belongs to the anaerobic coproporphyrinogen-III oxidase family.,

研究领域 (Research Area)

图片 (Image Data)




Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using RSAD1 Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from HeLa cells, using RSAD1 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: RSAD1 Rabbit Polyclonal Antibody
产品货号: APRab17400



Western Blot analysis of various cells using RSAD1 Polyclonal Antibody

注意事项 (Note)

For research use only .