

产品名称: S-100A7L2 Rabbit Polyclonal Antibody
产品货号: AP Rab17473



产品概述 (Summary)

产品名称 (Production Name)	S-100A7L2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	S100A7L2
别名 (Alternative Names)	S100A7L2; S100A7B; Protein S100-A7-like 2; S100 calcium-binding protein A7-like 2
基因 ID (Gene ID)	645922.0
蛋白 ID (SwissProt ID)	Q5SY68.The antiserum was produced against synthesized peptide derived from human S100A7L2. AA range:38-87

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	11kDa

产品名称: S-100A7L2 Rabbit Polyclonal Antibody
产品货号: APRab17473




研究背景 (Background)


similarity:Belongs to the S-100 family.,similarity:Contains 2 EF-hand domains.,similarity:Belongs to the S-100 family.,similarity:Contains 2 EF-hand domains.,

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of lysates from HUVEC cells, using S100A7L2 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HepG2 cells using S100A7L2 antibody.



产品名称: S-100A7L2 Rabbit Polyclonal Antibody
产品货号: APRab17473



Western Blot analysis of HuvEc cells using S-100A7L2 Polyclonal Antibody

注意事项 (Note)

For research use only .