

产品名称: S1A7A Rabbit Polyclonal Antibody
产品货号: APRab17492



产品概述 (Summary)

产品名称 (Production Name)	S1A7A Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	S100A7A S100A15 S100A7L1
别名 (Alternative Names)	
基因 ID (Gene ID)	338324.0
蛋白 ID (SwissProt ID)	Q86SG5.Synthesized peptide derived from human protein . at AA range: 30-110

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	11kDa

研究背景 (Background)

function:May be involved in epidermal differentiation and inflammation and might therefore be important for the


产品名称: S1A7A Rabbit Polyclonal Antibody
产品货号: APRab17492



pathogenesis of psoriasis and other diseases.,similarity:Belongs to the S-100 family.,similarity:Contains 2 EF-hand domains.,tissue specificity:Overexpressed in psoriasis.,function:May be involved in epidermal differentiation and inflammation and might therefore be important for the pathogenesis of psoriasis and other diseases.,similarity:Belongs to the S-100 family.,similarity:Contains 2 EF-hand domains.,tissue specificity:Overexpressed in psoriasis.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4°over night

注意事项 (Note)

For research use only .