

产品名称: SAA4 Rabbit Polyclonal Antibody
产品货号: AP Rab17568



产品概述 (Summary)

产品名称 (Production Name)	SAA4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SAA4
别名 (Alternative Names)	SAA4; CSAA; Serum amyloid A-4 protein; Constitutively expressed serum amyloid A protein; C-SAA
基因 ID (Gene ID)	6291.0
蛋白 ID (SwissProt ID)	P35542. The antiserum was produced against synthesized peptide derived from human SAA4. AA range: 81-130

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	

产品名称: SAA4 Rabbit Polyclonal Antibody
产品货号: APRab17568




研究背景 (Background)

function:Major acute phase reactant. Apolipoprotein of the HDL complex.,induction:Constitutively expressed.,similarity:Belongs to the SAA family.,tissue specificity:Expressed by the liver; secreted in plasma.,function:Major acute phase reactant. Apolipoprotein of the HDL complex.,induction:Constitutively expressed.,similarity:Belongs to the SAA family.,tissue specificity:Expressed by the liver; secreted in plasma.,

研究领域 (Research Area)

Cardiovascular; Blood; Acute Phase Reactants; Lipids / Lipoproteins; Lipoproteins/Apolipoproteins; Atherosclerosis; Lipoprotein metabolism; Apolipoproteins; Metabolism; Pathways and Processes; Metabolic signaling pathways; Lipid and lipoprotein metabolism; Cancer; Heart disease

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma, using SAA4 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .