

产品名称: SCAF1 Rabbit Polyclonal Antibody
产品货号: AP Rab17632



产品概述 (Summary)

产品名称 (Production Name)	SCAF1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SCAF1 SCAF1; SFRS19; SRA1; Splicing factor; arginine/serine-rich 19; SR-related and
别名 (Alternative Names)	CTD-associated factor 1; SR-related-CTD-associated factor; SCAF; Serine arginine-rich pre-mRNA splicing factor SR-A1; SR-A1
基因 ID (Gene ID)	58506.0
蛋白 ID (SwissProt ID)	Q9H7N4.The antiserum was produced against synthesized peptide derived from human SFRS19. AA range:41-90

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	140kDa

产品名称: SCAF1 Rabbit Polyclonal Antibody
产品货号: AP Rab17632





研究背景 (Background)

function: May function in pre-mRNA splicing., induction: Up-regulated by estrogens, androgens and glucocorticoids., similarity: Belongs to the splicing factor SR family., subunit: Interacts with POLR2A., tissue specificity: Ubiquitous. Highly expressed in fetal brain and liver, poorly expressed in salivary gland, heart, skin and ovary. Expressed in colorectal carcinomas and ovarian cancers. Overexpressed in colorectal carcinomas as compared to normal colonic mucosa., function: May function in pre-mRNA splicing., induction: Up-regulated by estrogens, androgens and glucocorticoids., similarity: Belongs to the splicing factor SR family., subunit: Interacts with POLR2A., tissue specificity: Ubiquitous. Highly expressed in fetal brain and liver, poorly expressed in salivary gland, heart, skin and ovary. Expressed in colorectal carcinomas and ovarian cancers. Overexpressed in colorectal carcinomas as compared to normal colonic mucosa.,

研究领域 (Research Area)

Epigenetics and Nuclear Signaling; DNA / RNA; RNA Processing Splicing

图片 (Image Data)



Western Blot analysis of various cells using SCAF1 Polyclonal Antibody diluted at 1: 2000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA) .

产品名称: SCAF1 Rabbit Polyclonal Antibody
产品货号: APRab17632



注意事项 (Note)

For research use only .