

产品名称: SCCA1/2 Rabbit Polyclonal Antibody
产品货号: APRab17639



产品概述 (Summary)

产品名称 (Production Name)	SCCA1/2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SERPINB3/SERPINB4
别名 (Alternative Names)	SERPINB3; SCCA; SCCA1; Serpin B3; Protein T4-A; Squamous cell carcinoma antigen 1; SCCA-1
基因 ID (Gene ID)	6317.0
蛋白 ID (SwissProt ID)	P29508.The antiserum was produced against synthesized peptide derived from the Internal region of human SERPINB3/4. AA range:111-160

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	45kDa

产品名称: SCCA1/2 Rabbit Polyclonal Antibody
产品货号: APRab17639




研究背景 (Background)


developmental stage:Its expression is closely related to cellular differentiation in both normal and malignant squamous cells.,function:May act as a protease inhibitor to modulate the host immune response against tumor cells.,similarity:Belongs to the serpin family.,similarity:Belongs to the serpin family. Ov-serpin subfamily.,subcellular location:Seems to also be secreted in plasma by cancerous cells but at a low level.,tissue specificity:Squamous cells. Expressed in some hepatocellular carcinoma (at protein level),,developmental stage:Its expression is closely related to cellular differentiation in both normal and malignant squamous cells.,function:May act as a protease inhibitor to modulate the host immune response against tumor cells.,similarity:Belongs to the serpin family.,similarity:Belongs to the serpin family. Ov-serpin subfamily.,subcellular location:Seems to also be secreted in plasma by cancerous cells but at a low level.,tissue specificity:Squamous cells. Expressed in some hepatocellular carcinoma (at protein level),,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysate from Jurkat cells, using SERPINB3/4 Antibody.




Western Blot analysis of Jurkat, 293 cells using SCCA1/2 Polyclonal Antibody.. Secondary antibody was diluted at 1:20000


产品名称: SCCA1/2 Rabbit Polyclonal Antibody
产品货号: APRab17639




Immunohistochemical analysis of paraffin-embedded human-liver-cancer, antibody was diluted at 1:100



Immunohistochemical analysis of paraffin-embedded human-liver-cancer, antibody was diluted at 1:100



Immunohistochemical analysis of paraffin-embedded human-lung, antibody was diluted at 1:100



产品名称: SCCA1/2 Rabbit Polyclonal Antibody
产品货号: APRab17639



Immunohistochemical analysis of paraffin-embedded human-lung, antibody was diluted at 1:100

注意事项 (Note)

For research use only .