

产品名称: SCCA2 Rabbit Polyclonal Antibody
产品货号: APRab17640



产品概述 (Summary)

产品名称 (Production Name)	SCCA2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SERPINB4
别名 (Alternative Names)	SERPINB4; PI11; SCCA2; Serpin B4; Leupin; Peptidase inhibitor 11; PI-11; Squamous cell carcinoma antigen 2; SCCA-2
基因 ID (Gene ID)	6318.0
蛋白 ID (SwissProt ID)	P48594.The antiserum was produced against synthesized peptide derived from human SERPINB4. AA range:261-310

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	45kDa

产品名称: SCCA2 Rabbit Polyclonal Antibody
产品货号: AP Rab17640




研究背景 (Background)

function: May act as a protease inhibitor to modulate the host immune response against tumor cells., similarity: Belongs to the serpin family., similarity: Belongs to the serpin family. Ov-serpin subfamily., subcellular location: Seems to also be secreted in plasma by cancerous cells but at a low level., tissue specificity: Squamous cells., function: May act as a protease inhibitor to modulate the host immune response against tumor cells., similarity: Belongs to the serpin family., similarity: Belongs to the serpin family. Ov-serpin subfamily., subcellular location: Seems to also be secreted in plasma by cancerous cells but at a low level., tissue specificity: Squamous cells.,


研究领域 (Research Area)

Immunology; Immune System Diseases; Allergens; Cell Biology; Proteolysis / Ubiquitin; Protease inhibitors; Serine protease inhibitors; SERPINs; Cancer; Tumor biomarkers; Enzymes

图片 (Image Data)




Western blot analysis of SERPINB4 Antibody. The lane on the right is blocked with the SERPINB4 peptide.



Western blot analysis of the lysates from COLO205 cells using SERPINB4 antibody.

产品名称: SCCA2 Rabbit Polyclonal Antibody
产品货号: APRab17640



Western Blot analysis of various cells using SCCA2 Polyclonal Antibody

注意事项 (Note)

For research use only .