

产品名称: SCFD1 Rabbit Polyclonal Antibody
产品货号: APRab17643



产品概述 (Summary)

产品名称 (Production Name)	SCFD1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SCFD1
别名 (Alternative Names)	SCFD1; C14orf163; KIAA0917; STXBP1L2; FKSG23; Sec1 family domain-containing protein 1; SLY1 homolog; Sly1p; Syntaxin-binding protein 1-like 2
基因 ID (Gene ID)	23256.0
蛋白 ID (SwissProt ID)	Q8WVM8.The antiserum was produced against synthesized peptide derived from human SCFD1. AA range:533-582

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	

产品名称: SCFD1 Rabbit Polyclonal Antibody
产品货号: APRab17643




研究背景 (Background)

function:Involved in vesicular transport between the endoplasmic reticulum and the Golgi.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the STXBP/unc-18/SEC1 family.,subunit:Binds STX5A.,function:Involved in vesicular transport between the endoplasmic reticulum and the Golgi.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the STXBP/unc-18/SEC1 family.,subunit:Binds STX5A.,


研究领域 (Research Area)

Signal Transduction; Protein Trafficking; Vesicle Transport; Regulation; ER Proteins; Golgi Proteins

图片 (Image Data)




Western blot analysis of SCFD1 Antibody. The lane on the right is blocked with the SCFD1 peptide.



Western blot analysis of the lysates from 293 cells using SCFD1 antibody.

产品名称: SCFD1 Rabbit Polyclonal Antibody
产品货号: APRab17643



Western Blot analysis of various cells using SCFD1 Polyclonal Antibody

注意事项 (Note)

For research use only .