

产品名称: Secretin Rabbit Polyclonal Antibody
产品货号: APRab17698



产品概述 (Summary)

产品名称 (Production Name)	Secretin Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SCT
别名 (Alternative Names)	Secretin
基因 ID (Gene ID)	6343.0
蛋白 ID (SwissProt ID)	P09683.Synthesized peptide derived from human Secretin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	13kDa

研究背景 (Background)

This gene encodes a member of the glucagon family of peptides. The encoded preproprotein is secreted by endocrine S

产品名称: Secretin Rabbit Polyclonal Antibody
产品货号: APRab17698




cells in the proximal small intestinal mucosa as a prohormone, then proteolytically processed to generate the mature peptide hormone. The release of this active peptide hormone is stimulated by either fatty acids or acidic pH in the duodenum. This hormone stimulates the secretion of bile and bicarbonate in the duodenum, pancreatic and biliary ducts. [provided by RefSeq, Feb 2016],function:Stimulates formation of NaHCO_3 -rich pancreatic juice and secretion of NaHCO_3 -rich bile and inhibits HCl production by the stomach.,online information:Secretin entry,similarity:Belongs to the glucagon family.,

研究领域 (Research Area)

Signal Transduction; Growth Factors/Hormones; Hormones; Tags & Cell Markers; Cell Type Markers; Other Cell Types

图片 (Image Data)



Western blot analysis of various lysates, primary antibody was diluted at 1:1000, 4° over night, secondary antibody was diluted at 1:10000, 37° 1 hour.

注意事项 (Note)

For research use only .