

产品名称: Septin 14 Rabbit Polyclonal Antibody
产品货号: APRab17739



产品概述 (Summary)

产品名称 (Production Name)	Septin 14 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SEPT14
别名 (Alternative Names)	SEPT14; Septin-14
基因 ID (Gene ID)	346288.0
蛋白 ID (SwissProt ID)	Q6ZU15.The antiserum was produced against synthesized peptide derived from human SEPT14. AA range:341-390

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	50kDa

研究背景 (Background)

产品名称: Septin 14 Rabbit Polyclonal Antibody
产品货号: APRab17739




SEPT14 is a member of the highly conserved septin family of GTP-binding cytoskeletal proteins implicated in membrane transport, apoptosis, cell polarity, cell cycle regulation, cytokinesis, and other cellular functions (Peterson et al., 2007 [PubMed 17922164]).[supplied by OMIM, Jul 2008],function:Involved in cytokinesis .,similarity:Belongs to the septin family.,


研究领域 (Research Area)

Cell Biology; Cell Cycle; Cell Division; Cytokinesis


图片 (Image Data)



Western blot analysis of SEPT14 Antibody. The lane on the right is blocked with the SEPT14 peptide.



Western blot analysis of the lysates from Jurkat cells using SEPT14 antibody.



Western Blot analysis of 3T3 cells using Septin 14 Polyclonal Antibody

产品名称: Septin 14 Rabbit Polyclonal Antibody
产品货号: APRab17739



注意事项 (Note)

For research use only .