

产品名称: Septin 6 Rabbit Polyclonal Antibody
产品货号: AP Rab17743



产品概述 (Summary)

产品名称 (Production Name)	Septin 6 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SEPT6
别名 (Alternative Names)	SEPT6; KIAA0128; SEP2; Septin-6
基因 ID (Gene ID)	23157.0
蛋白 ID (SwissProt ID)	Q14141.The antiserum was produced against synthesized peptide derived from human SEPT6. AA range:223-272

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	50kDa

研究背景 (Background)

产品名称: Septin 6 Rabbit Polyclonal Antibody
产品货号: AP Rab17743




septin 6(SEPT6) Homo sapiens This gene is a member of the septin family of GTPases. Members of this family are required for cytokinesis. One version of pediatric acute myeloid leukemia is the result of a reciprocal translocation between chromosomes 11 and X, with the breakpoint associated with the genes encoding the mixed-lineage leukemia and septin 2 proteins. This gene encodes four transcript variants encoding three distinct isoforms. An additional transcript variant has been identified, but its biological validity has not been determined. [provided by RefSeq, Jul 2008],function:Involved in cytokinesis .,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the septin family.,subunit:Heterohexamer composed of two heterotrimers containing one copy each of SEPT2, SEPT6 and SEPT7. The asymmetrical heterotrimers associate head-to-head to form a hexameric unit that assembles into filaments. Interacts with SEPT12. Interaction with SEPT12 alters filament structure.,


研究领域 (Research Area)

Cell Biology; Cell Cycle; Cell Division; Cytokinesis; Microbiology; Interspecies Interaction; Host Virus Interaction

图片 (Image Data)




Western blot analysis of the lysates from HeLa cells using SEPT6 antibody.



Western Blot analysis of JK cells using Septin 6 Polyclonal Antibody diluted at 1: 1000

产品名称: Septin 6 Rabbit Polyclonal Antibody
产品货号: APRab17743



Western blot analysis of various lysis using Septin 6 Polyclonal Antibody diluted at 1: 1000. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .