

产品名称: Serinc1 Rabbit Polyclonal Antibody
产品货号: APRab17752



产品概述 (Summary)

产品名称 (Production Name)	Serinc1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SERINC1
别名 (Alternative Names)	SERINC1; KIAA1253; TDE1L; TDE2; Serine incorporator 1; Tumor differentially expressed protein 1-like; Tumor differentially expressed protein 2
基因 ID (Gene ID)	57515.0
蛋白 ID (SwissProt ID)	Q9NRX5.The antiserum was produced against synthesized peptide derived from human SERC1. AA range:381-430

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	48kDa

产品名称: Serinc1 Rabbit Polyclonal Antibody
产品货号: APRab17752




研究背景 (Background)

function:Enhances the incorporation of serine into phosphatidylserine and sphingolipids.,similarity:Belongs to the TDE1 family.,subunit:Interacts with SPTLC1.,function:Enhances the incorporation of serine into phosphatidylserine and sphingolipids.,similarity:Belongs to the TDE1 family.,subunit:Interacts with SPTLC1.,


研究领域 (Research Area)

Lipid metabolism; Metabolism; Pathways and Processes; Metabolic signaling pathways; Amino acid metabolism; Amino Acids; Signal Transduction

图片 (Image Data)



Western blot analysis of lysates from HepG2 cells, using SERC1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using Serinc1 Polyclonal Antibody

注意事项 (Note)

For research use only .